Bite registration material, silicone based, addition curing fast setting – extreme high final hardness


Fig. 1

1. Mixing and Dosing

Extrude material from the cartridge with Automix2-system gun (fig. 1). Place cartridge into mixing gun and remove closure cap by rotation. Eject a small amount of material (bleeding) until equal flow is being extruded from both orifices (fig. 2). Insert mixing cannula (green) into guide grooves on cartridge and cannula and rotate anti-clockwise to a stop (fig. 3). Store cartridge with mixing cannula attached for subsequent use.


Fig. 2

2. Instruction for use

Extrude material from cartridge with mixing cannula or with attached contouring tip onto teeth (fig. 4). At terminal occlusion application from vestibular is possible. Let registration harden in terminal occlusion and remove from the mouth after 1 min.

Fig. 3

3. Disinfection

After removal from the mouth rinse registrate under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.

Important working hints

- Remove any retraction astringents or mouth wash solutions prior to the registration – could affect setting reaction.
- Do not intermix with condensation curing silicones.
- Avoid staining of clothing. Hardened impression materials cannot be removed.
- Remove all residues of registration material from the patient's mouth.
- For reducing the pressure needed during injection, we recommend using the green mixing cannula.


Fig. 4

Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

For use by trained specialists.

Indications for use:

- Bite registration with crown and bridge impressions and for orthodontic indications, key for intraoral registration

Technical Data:

- DIN 13903
- **Mixed volume:** 50 ml (cartridge)
- **Mixing ratio:** 1:1
- **Colour code:** base: green catalyst: yellow
- **Mixing time:** omitted (Automix2-system)
- **Working time:** approx. 30 sec.*
- **Setting time:** approx. 1 min. 30 sec.*
- **Time in mouth:** approx. 1 min.
- **Hardness Shore D:** 45
- **Linear dimensional change:** < 0.2 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity
- **Storage:**

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

Ordering information:

greenbite apple
Standard packing **02641**
2 x 50 ml cartridges
12 mixing cannulas, green
12 contouring tips for mixing cannulas green

Ecopacking **02647**
8 x 50 ml cartridges
48 mixing cannulas, green
48 contouring tips for mixing cannulas green

Mixing cannulas **02705**
green, 48 pcs.

Contouring tips **02648**
for mixing cannulas green
12 pcs.

milkbite
Standard packing **03380**
2 x 50 ml cartridges
12 mixing cannulas, pink
12 contouring tips

greenbite colour
Standard packing **02207**
2 x 50 ml cartridges
12 mixing cannulas, green
12 contouring tips

glassbite
Standard packing **02571**
1 x 50 ml cartridge
4 mixing cannulas, green
4 contouring tips

pixelbite®
Standard packing **03103**
2 x 50 ml cartridges
12 mixing cannulas, pink
12 contouring tips

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