3M ESPE


- Ketac™ Cem radiopaque**
- Glass Ionomer Luting Cement**
- Рентгенконтрастный стеклоиономерный цемент для фиксации**
- Глас-иономерен цимент за залепване**
- Staklenoiomerni cement za definitivno cementiranje**
- Üvegionomer ragasztó cement**
- Szko-jonomerowy cement do osadzania**
- Ciment ionomer de sticlă pentru cimentări**
- Skloionomérný tmeliaci cement**
- Glasionomerni cement za stalno cementiranje**
- Skloionomerni tmelici cement**
- Cam ionomer luting sement**
- Klaasionomeerne kinnitustsemnt**
- Stikka jonomëra lutinga cements**
- Stiklo jonomerinis fiksavimo cementas**

Instructions for Use
Инструкции по применению
Указания за употреба
Upute za upotrebu
Kullanma Talimatları
Kasutusalutused
Instrukcja użycia
Instruțiuni de Utilizare

Made in Germany by
3M ESPE AG
Dental Products
D-82229 Seefeld – Germany

3M ESPE
Dental Products
St. Paul, MN 55144-1000

3M, ESPE, Alkaliner and Ketac are trademarks of 3M or 3M ESPE AG.
© 2003. All rights reserved.

4400763646401

ENGLISH

Product Description
Ketac™ Cem radiopaque, manufactured by 3M ESPE, is a fixing cement for mixing by hand.

For details on Alkaliner™, manufactured by 3M ESPE, please refer to corresponding Instructions for Use. The Instructions for Use on the product concerned must be kept available for the length of time it is employed.

Indication

- Cementation of inlays, onlays, crowns and bridges
- Luting cement for use on orthodontic strips
- Relinings

Preparation

- To provide for optimal adhesion, carefully clean enamel, dentin, and metal surfaces.
- Blow the cavity dry in only 2-3 short intervals with water and oil-free air or dab dry with cotton pellets. Do not over-dry!
 - The cavity should be just so dry that the surface has a matt shiny appearance.** Excessive drying with compressed air can negatively affect adhesion or result in post-operative sensitivity after seating the restoration.
- The metal surfaces should be dried well with compressed air.
- Avoid further contamination by saliva.

Pulp Protection
The glass ionomer cement must not be applied directly to dentin – e.g. on inlays – situated close to the pulp or to the exposed pulp. Prior to taking the impression, cover areas close to the pulp using a hard-setting calcium hydroxide preparation, e.g. Alkaliner.

Dosage

- Shake bottle to loosen powder.
- Screw off cap and replace by dispensing dropper.
- For cementing work, the recommended powder-to-liquid mixing ratio is 3.8:1 w/w, this corresponding to 1 level spoonful of powder to two drops of liquid.
- Skim off spoon on plastic segment. Do not compress the powder.
- Dose powder and liquid next to each other on a mixing pad or glass slab.
 - Keep bottle upright while dispensing.
 - Crystallized liquid must not be allowed to form on the dispensing dropper.
- For relinings more powder may be added according to the desired consistency.
- Carefully seal bottle after use.

If the mix is too thin (due to insufficient dosing of powder) this may produce symptoms of pulpitis.

Mixing
Working of Ketac Cem radiopaque should proceed at a room temperature of 20-25°C/68-77°F.

- Mix using a metal or plastic cement spatula.
- Add the powder to the liquid in one portion.
 - On principle a sufficient quantity of powder should be added to the liquid to produce viscid consistency.
 - The quantity of cement required for a crown should just about drip from the spatula.
- Continue to smooth out the paste until a homogeneous mixture is obtained.

Application

- During application, water and saliva must be prevented from entering the working area.
- Apply a thin coat of cement to the inside of the crown and stump.
- Afterwards proceed with restoration.
 - Avoid overfilling the crown. Applying hydrostatic pressure buildup in closely adjoining crowns may cause pulp discomfort.

Times
The following time scale applies at a room temperature of 23°C/73°F and 50% relative atmospheric humidity:

	min:sec
Mixing	1:00
Application incl. mixing	3:30
Setting from begin of mix	7:00

Higher temperatures shorten the working time, whereas lower temperatures will extend this period (e.g. mixing on a cooled glass slab). A higher amount of powder also brings about a shorter working time. Exceeding the working time causes loss of adhesion to enamel and dentine.

Removal of Excess

- Remove excess using a Heidemann spatula and/or probe 7 min. after beginning of mix.

Incompatibilities
In susceptible individuals, sensitization to the product cannot be excluded. If allergic reactions are observed, use of the product should be discontinued and the product completely removed. In rare cases, pulp complications could occur, especially if the recommendations for processing are disregarded.

Storage and Shelf Life
Keep the powder in a dry place away from damp. Do not store the product above 25°C/77°F. Do not use after the expiry date.

Customer Information
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty
3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user’s application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE’s sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability
Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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