ENGLISH

General Information

RelyX™ Veneer Cement, manufactured by 3M ESPE, is a radiopaque, color stable, light-cured resin cement indicated for bonding veneers fabricated of porcelain or composite. It is used in conjunction with RelyX™ Try-In Pastes, Adper™ Single Bond 2 Adhesive or Adper™ Scotchbond™ Multi-Purpose Adhesive and RelyX™ Ceramic Primer, manufactured by 3M ESPE. The RelyX Veneer cement contains a resin system consisting of a TEGDMA/BisGMA blend and a zirconia/silica filler with a particle size range of 0.2–3.0 microns and volume percentage of approximately 47%.

Precautions

RelyX Veneer cement contains acrylate resins. Avoid use of this product on patients with known acrylate allergies. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured resins. **Use of protective gloves and a no-tooth technique is recommended.**

If skin contact occurs, wash skin with soap and water. Acrylates may penetrate commonly-used gloves. If cement contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If accidental contact with eyes or prolonged contact with oral soft tissue occurs, flush with large amounts of water. If irritation persists, consult a physician.

Instructions for Use

RelyX Veneer cement is designed to be used in conjunction with Adper Single Bond 2 Adhesive or Adper Scotchbond Multi-Purpose adhesive and RelyX Try-In pastes. It is recommended that RelyX Veneer cement be used only with 3M ESPE adhesive systems where veneer bonding is indicated. Please read the instructions for use of the 3M ESPE adhesive system you use for additional safety and dispensing information. **Note: 3M ESPE disclaims liability for claims arising from the use of RelyX Veneer cement with any adhesive other than Adper Single Bond 2 Adhesive and Adper Scotchbond Multi-Purpose Adhesive.**

1. To prepare the veneer, the internal bonding surface of a ceramic veneer surface should be microetched with a hydrofluoric acid etchant. The dental laboratory typically does this. The internal surface of a composite veneer should be roughened by sandblasting or air abrasion.

2. Remove the provisional restoration, if present, from the tooth. Clean tooth with pumice; rinse thoroughly and lightly air dry.

3. Check the fit and esthetics of the restoration by trying in the veneer with water or the Translucent shade of RelyX Try-In paste. If the fit and esthetics are acceptable, proceed directly to step 5 using the Translucent shade of RelyX Veneer cement.

4. If color adjustment is needed, select the appropriate shade of RelyX Try-In paste. Apply RelyX Try-In paste directly onto the internal surface of the veneer. Gently seat the veneer onto the tooth preparation. With multiple veneers, start at the midline and seat veneers sequentially working back. The shade of RelyX Try-In paste is formulated to precisely match the corresponding cured shade of the cement. Once the fit and esthetics are verified, remove the veneers and thoroughly rinse the RelyX Try-In paste from the tooth and restoration with a water spray and then dry.

5. The internal surface of the veneer should be thoroughly cleaned prior to applying the silane to remove any contaminants. The surface can be cleaned by applying Scotchbond™ Etchant, manufactured by 3M ESPE (35% phosphoric acid) for 15 seconds followed by rinsing and drying. Dispense RelyX Ceramic Primer (silane) into a mixing well. Brush a layer of the silane onto the internal surface of the veneer. Lightly air-dry the surface for 5 seconds to evaporate the solvent.

6. Clean the prepared teeth using flour of pumice slurry. Rinse and lightly dry. Isolate the area to prevent contamination. Place matrix strips between teeth to prevent unwanted bonding to adjacent teeth.

7. Etch the tooth preparation. Apply Scotchbond etchant to both dentin and enamel for 15 seconds. Rinse for 10 seconds. **Blot excess water leaving tooth moist. If the tooth is accidentally dried, rewet the bonding surface with water for 15 seconds and blot excess leaving tooth moist.** Alternatively, a HEMA based desensitizing agent can be applied to the etched dentin according to the manufacturer's instructions.

8. Apply the adhesive. **For Adper Single Bond 2 Adhesive:** Immediately after blotting, apply 2–3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. **Do not light cure.** Apply 1 coat of Adper Single Bond 2 Adhesive to the silane treated veneer. Dry gently for 5 seconds. **Do not light cure. For Adper Scotchbond Multi-Purpose adhesive:** Apply Adper Scotchbond Multi-Purpose Primer to etched enamel and dentin. Dry gently for 5 seconds. Apply Adper Scotchbond Multi-Purpose Adhesive to the primed enamel and dentin. **Do not light cure.** Apply 1 coat Adper Scotchbond Multi-Purpose Adhesive to the silane treated veneer. **Do not light cure.**

9. Apply a thin layer of the selected shade of RelyX Veneer cement directly from the syringe onto the bonding surface of the veneer. **Protect the cement from direct operatory light exposure.**

10. Seat the veneer in place using gentle pressure. Excess cement will extrude out around the margins. Spot cure the veneer into place on the facial surface away from the margins using a small diameter light guide for 20 seconds to secure it in place.

11. Remove the excess cement from the margins using a blunt instrument or dry brush. Alternatively, the excess cement can be "back-cured" for 3 seconds **with a standard halogen light** to allow the excess to gel and allow for easier cleanup.

12. Light cure the labial, lingual, interproximal and occlusal surfaces for 30 seconds each. The A5 Opaque/Dark shade will require a 40 second exposure for each surface. Curing recommendations will allow for cement cure depths not less than 0.5mm for the opaque shades and not less than 1.5 mm for the non-opaque shades. **Note: The cure times are indicated for a standard halogen curing light with a minimum light output of 400 mW/cm².**

13. Remove matrix strips. Finish marginal areas with Sof-Lex™ Finishing and Polishing Strips and Discs, manufactured by 3M ESPE. Check and adjust occlusion as necessary.

Storage and Use

1. RelyX Veneer cement is designed to be used at room temperature of 2–27° C or 35–80° F.
2. Shelf life in a refrigerator or at room temperature is 24 months. See outer package for expiration date. Storage in refrigerator ensures longest possible shelf life. Allow cement to come to room temperature before use.
3. Do not store in proximity to eugenol-containing products.

4. Do not expose the cement to elevated temperatures or intense light.

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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