**Precision impression material, C-silicone based,  
medium consistency – medium bodied**

Fig. 1


Fig. 2

## 1. Preparation of the impression

Depending on the impression technique select an appropriate tray. For optimal adhesion of the impression to the tray apply a thin layer of **sili Adhesive Spray** or **sili Adhesive** liquid onto the tray and let dry for 60 sec. (Fig. 1).

## 2. Mixing

Extrude **silone®** and **cat<sup>p</sup> universal FUTUR** from tubes in equal length onto mixing pad (Fig. 2). With dosage of liquid catalyst 1 cm **silone®** to 1 drop of catalyst. Mixing 30 sec. with spatula until homogeneous mix and uniform colour is achieved.

## 3. Disinfection

After removal rinse impression under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.

## 4. Model casting

Casting between 30 minutes and 24 hours after removal of the impression. Recommended materials: Dental plasters of class III.

## 5. Galvanisation

Impressions can be electroplated with commercially available copper and silver plating baths.

## 6. Cleaning of trays

Remove set material with a blunt instrument. **sili Spray** or **sili liquid** film can be easily removed with **sili Spray Solvent** (free of FCHC). Work in well ventilated rooms. Clean and disinfect as usual.

## Important working hints


- Increased temperatures accelerate, decreased temperatures retard the setting time.
- Do not intermix with addition curing silicones.
- Catalyst bottles or tubes should be closed tightly immediately after use.
- Don't leave any residual material in the patient's mouth.
- Avoid staining of clothing. Cured impression materials cannot be removed.
- Avoid contact of catalyst with eyes, may cause irritation. If necessary flush eyes with plenty of water and seek immediate medical care.

Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

For use by trained specialists.

## silone® / cat<sup>p</sup> universal FUTUR:

## Indications for use:

- Universal material for dentulous impression
- Partial and copping and inlay impressions

## Technical Data:

DIN EN ISO 4823 – Type 2

### ■ Mixing ratio: silone® / cat<sup>p</sup> universal FUTUR

1 cm silone® = 1.00 g /  
1 cm cat<sup>p</sup> = 0.05 g

**silone® / sta-seal cat<sup>f</sup>**  
1 cm silone® = 1.10 g /  
1 drop cat<sup>f</sup> = 0.02 g

### ■ Colour code: silone®: red-brown cat<sup>p</sup> universal FUTUR: blue

**sta-seal cat<sup>f</sup>:**  
light blue

### ■ Mixing time: approx. 30 sec.

### ■ Working time: approx. 1 min. 15 sec.\*

### ■ Setting time: approx. 3 min. 45 sec.\*

### ■ Time in mouth: approx. 2 min. 30 sec.

### ■ Strain in compression: 9.50 %

### ■ Recovery from deformation: 98.70 %

### ■ Linear dimensional change: 0.45 %

### ■ Application: At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity

### ■ Storage:

15°C  
59°F

25°C  
77°F

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

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kneadable,  
very high  
consistency

## silaplast FUTUR

high final hardness  
normal setting

Preliminary impression  
material for corrections of  
preparations, tray material,  
overimpressions, wax-bite

Wash  
material  
light bodied

low viscosity  
low con-  
sistency

## silasoft® N

normal setting

Correction impression  
material for the double  
impression and double  
mix technique, relinings

Wash  
material  
light bodied

low viscosity  
low con-  
sistency

## silasoft® S

normal setting

Correction impression  
material for the double  
impression and double  
mix technique, relinings

Functional  
heavy bodied

high con-  
sistency and  
plasticity

## sta-seal f

medium final hardness  
plastic phase  
extended setting time

Functional impressions,  
functional border moulding

Catalyst  
paste

Catalyst  
liquid

## cat<sup>p</sup> / cat<sup>f</sup> FUTUR

cat<sup>p</sup> universal paste,  
cat<sup>f</sup> silaplast liquid,  
sta-seal cat<sup>f</sup>, liquid

Special  
instrument

Accessory

## Deta-Cut

## sili adhesive

## sili solvent

Sterilizable, curved, for  
cutting the preliminary  
impression

adhesive for C-silicones

removal of the adhesive  
and cleaning of trays

All DETAX condensation-curing silicones are portioned with dispensing aids (tray, measuring bowl, measuring scale etc.) and kneaded by hand or mixed homogeneously on a mixing pad.

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