**Precision impression material, silicone based, condensation curing, medium consistency – medium bodied**

Fig. 1

Fig. 2

## 1. Preparation of the impression

Depending on the impression technique select an appropriate tray. For optimal adhesion of the impression to the tray apply a thin layer of **sili Adhesive** spray or **sili Adhesive** fluid onto the tray and let dry for 60 sec. (Fig. 1).

## 2. Mixing

Extrude **sta-seal** and **cat<sup>®</sup> universal FUTUR** from tubes in equal length and diameter onto mixing pad (Fig. 2). With dosage of liquid catalyst 1 cm **sta-seal** to 1 drop of catalyst. Mixing 30 sec. with spatula until homogeneous mix and uniform colour is achieved.

## 3. Disinfection

After removal rinse impression under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.

## 4. Model casting

Casting between 2 hours and 24 hours after removal of the impression. Recommended materials: Dental plasters of class III.

## 5. Galvanisation

Impressions can be electroplated with commercially available copper and silver plating baths.

## 6. Cleaning of trays

Remove set material with a blunt instrument. **sili Adhesive** film can be easily removed with **sili Solvent** spray (free of FCHC). Work in well ventilated rooms. Clean and disinfect as usual.

## Important working hints


- Overdosing of catalyst or increased temperatures accelerate, underdosing or decreased temperatures retard the setting time.
- Catalyst bottles or tubes should be closed tightly immediately after use.
- Avoid staining of clothing. Hardened impression materials cannot be removed.
- Avoid contact of catalyst with eyes, may cause irritation. If necessary flush eyes with plenty of water and seek medical care immediately.
- Don't leave any residual material in the patient's mouth.

Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to get tested the allergenic potential before the application of the material.

For use by trained specialist staff.

## sta-seal / cat<sup>®</sup> universal FUTUR:

## Indications for use:

- Universal impression material for edentulous impressions
- Mucostatic and relined impression

## Technical Data:

DIN EN ISO 4823 – Type 2

### Mixing ratio:

**sta-seal / cat<sup>®</sup> universal FUTUR**

1 cm sta-seal = 1.10 g / 1 cm cat<sup>®</sup> = 0.05 g

**sta-seal / sta-seal cat<sup>®</sup>**  
1 cm sta-seal = 1.10 g / 1 drop cat<sup>®</sup> = 0.02 g

### Colour code:

**sta-seal:** lime green  
**cat<sup>®</sup> universal FUTUR:** blue

**sta-seal cat<sup>®</sup>:** light blue

### Mixing time:

approx. 30 sec.

### Working time:

1 min. 15 sec.\*

### Setting time:

approx. 3 min. 45 sec.\*

### Time in mouth:

approx. 2 min. 30 sec.

### Strain in compression:

14.50 %

### Recovery from deformation:

99.30 %

### Linear dimensional change:

0.45 %

### Application:

At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity

### Storage:

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

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medium bodied  
medium consistency

## silone®

normal setting  
final hardness approx.  
45 Shore A

Universal material for all dentulous impressions, single phase impression technique, for partial-, coppering and inlay impressions

medium bodied  
medium consistency

## sta-seal

normal setting  
final hardness approx.  
35 Shore A

Precision impression material for edentulous impressions, atrophic conditions, flabby ridges, mucostatic and relined impressions

heavy bodied  
high consistency & plasticity

## sta-seal f

ideal final hardness  
extended plastic phase  
extended setting time

Special impression material for functional impressions with extended setting, functional border moulding, functional impression in one session with si-plast trays

catalyst paste  
catalyst liquid

## cat<sup>p</sup> / cat<sup>f</sup> FUTUR

inert system  
non-irritating

cat<sup>p</sup> universal paste,  
silasoft® cat<sup>f</sup>, liquid,  
sta-seal cat<sup>f</sup>, liquid

special instrument  
accessory

## Deta-Cut

## sili adhesive

## sili solvent

Sterilizable, curved, for cutting the preliminary impression

adhesive for C-silicones

removal of the adhesive and cleaning of trays

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