tempolink® single/multi

GB
Instructions for use

Temporary luting cement for crowns and bridges,
eugenol-free – polycarboxylate-based


Fig. 1

1. Cleaning the preparation

Before cementing, clean the preparations with water. If the preparation has been disinfected beforehand, remove all traces of disinfectant from the surface. Dry the preparation area (Fig. 1).


Fig. 2

2. Mixing and applying

Dispensing is performed with the mini-mix system. Insert plunger into the body of the cartridge. Remove the cap from the cartridge by turning. Before fitting the mixing cannula, dispense a small quantity of material until material emerges from the two openings uniformly (Fig. 2). When attaching the mixing cannula note the guides on the cannula and cartridge. Lock in place by turning in the opposite direction (Fig. 3) – before each application dispense and discard a small quantity. Dispensing can now be performed to suit requirements. Dispense material with uniform pressure. After use, leave the mixing cannula on the cartridge until the next application.

Fig. 3

3. Inserting the temporary restoration

Apply a thin layer of **tempolink®** to the interior of the temporary restorations (Fig. 4). Directly after that, position the temporary restorations on the preparations by applying slight pressure (Fig. 5), then press them on firmly. Remove any excess material either immediately after insertion using a moist cotton swab or before reaching the setting time using a suitable instrument (scaler or dental floss, Fig. 6).

Fig. 4

4. Working time

Working time at room temperature:
tempolink® single 60 sec.
tempolink® multi 2 min. 30 sec.

Curing in the mouth at body temperature:
tempolink® single 60 - 90 sec.
tempolink® multi 3 min. - 3 min. 30 sec.
Setting outside the mouth can be delayed.


Fig. 5

Important working hints

- **tempolink®** is compatible with permanent cements.
- **tempolink®** does not have a delaying effect on the curing of composite cement or composite restorative material.
- Due to the acidity of one of the components, pain may occur in case of preparations with open dentinal tubules. Therefore, we recommend to apply **smartprotect®** or **smartprotect® soft** prior to the insertion of the restoration.
- The cement is not allowed to be applied directly in close proximity to the unprotected pulp.
- Cement which has set can be removed from instruments etc. using 10% sodium hydroxide solution.
- Manual removal and mixing of the components is not recommended.
- Do not allow individual components to come into contact with skin or mucous membrane.
- Using powder/monomer materials, wait the complete setting and roughen the provisionals on the inside.
- Do not leave any residues of material in the mouth.




Fig. 6

For dental use only.

tempolink® single

tempolink® multi

Applications:

- Temporary luting of crowns, bridges, inlays & onlays
- Temporary luting of permanent restorations
- **tempolink® single:** single and small, temporary restorations
- **tempolink® multi:** large temporary restorations

Technical data:

- **Mixed volume:** 5 ml (mini-mix)
- **Mixing ratio:** 1:1
- **Colour code:** Base: natural-opaque
Catalyst: transparent
- **Mixing time:** omitted (mini-mix)
- **Working time:**
single: approx. 60 sec.
multi: approx. 150 sec.
- **Setting time in the mouth:**
single: approx. 1 min. - 1 min. 30 sec.
multi: approx. 3 min. - 3 min. 30 sec.
- **Application:**
At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity
- **Storage:**
At 15-25 °C / 59-77 °F, < 50 % rel. humidity

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

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