VOCO Twinky Star

GB	CE ⁰⁴⁸²
in accordance with DIN EN ISO 4049	Instructions for use
DE	entspricht DIN EN ISO 4049
FR	correspond à DIN EN ISO 4049
ES	corresponde a DIN EN ISO 4049
PT	En conformidade com a DIN EN SO 4049
IT	In corrispondenza con DIN EN ISO 4049
RU	Используется по инструкции
PL	Wskazywanie sposobu użycia

Product description:
Twinky Star is a coloured, light-curing and radiopaque compomer filling system with glitter effect for cavities in deciduous teeth.

Twinky Star can be cured with halogen light/LED light (blue light) and is used with a dentine/enamel bond.

Indications:
- fillings in deciduous teeth

Use:
Preparation/Colour Selection
Clean teeth to be treated with fluoride-free cleaning paste. Choose one of the 8 colours indicated in the shade guide. Secure a dry working field.

Preparation of the cavity:
Principally, only a minimal preparation (adhesive filling technique) must be carried out in order to preserve healthy tooth substance. Follow special preparation forms for deciduous teeth!

Cleaning:
- Remove any residues in the cavity with a water jet.
- Avoid contamination with blood or saliva after cleaning.
- Remove any excess water with a faint airjet. Do not overdry the dentine!

The dentine surface should be moist but not too wet (wet bonding).

Pulp protection:
It is recommended to apply a calcium hydroxide lining in areas close to the pulp.

Bonding agents:
For optimal adhesion, a conventional dentine/enamel bond should be used according to manufacturer’s instructions.

Application of Twinky Star:
Fillings over 2 mm in depth should be applied and cured in layers. Each layer must be polymerised for approx. 40 s. To light-cure this material, conventional polymerisation devices are suited. The light output should be a minimum of 500 mW/cm² on halogen polymerisation devices and 300 mW/cm² on LED devices. Use transparent strips and crowns or light-wedges. Ensure a tight marginal seal. The tip of the light curing device should be held as close as possible to the filling surface. If this distance is more than 5 mm, curing depth may be reduced or curing limited to the light cone.
Insufficient curing may lead to discolouration or pulpitis-like complaints. Finishing and polishing (with cooling) can be carried out immediately after removal of the shaping aids (e.g. with ultra-fine diamonds, polishing disks). The filling margin or preferably the entire tooth should be fluoridated as a finishing step.

Warnings and precautionary measures:
- **Twinky Star** contains Bis-GMA, diurethane dimethacrylate, TEGDMA, BHT. **Twinky Star** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- **Twinky Star** must not be used for direct or indirect pulp capping or occlusion-bearing, permanent fillings of permanent teeth.
- Preparations containing eugenol or clove oil disturb polymerisation of synthetic restoratives. Do not use zinc oxide eugenol cements or other materials containing eugenol when applying **Twinky Star**.
- **Twinky Star** is a compomer-based restorative with special colour effects. The glitter effect is caused by an insoluble mineral. **Twinky Star** does not contain any metals.
All colours comply with food and cosmetics regulations and are toxicologically safe. They are added in very low concentrations and are insoluble. After polymerisation, the colours are firmly bound in the synthetic matrix. **Twinky Star** is safe to use.

Storage:
Store at temperatures of 4 °C to 23 °C.
Allow **Twinky Star** to reach room temperature before use.
Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 AA1680 E1 0712 99 © by VOCO